




The Olympus Xenon E180 light source offers economical access to superior Xenon technology.

- Bright image illumination by high intensive 180 W Xenon lamp
- Excellent reliability due to lamp life counter
- Naturally color recognition as a result of the daylight resembling Xenon light


## Light Source Xenon E180

WA58675A Light source "Xenon E180",  
xenon,  
180 W



## Spare Lamp

WA58677A Lamp,  
xenon, 180 W,  
for WA58675A,  
2 pcs.



## Technical Data

### Power supply

Voltage ..... 100 – 240 V

Frequency ..... 50/60 Hz

Power consumption ..... 220 W

### Size

Dimensions ..... 275 (W) x 110 (H) x 327 (D) mm

Weight ..... 6 kg

### Classification

Protection against electric shock ..... class 1  
type BF

### Lamp

Type ..... Xenon

Power ..... 180 W

Lamp life ..... at least 500 hrs.<sup>1</sup>

<sup>1</sup> in intermittent use test cycle (1 hr. on, 1 hr. off)